


	<b>Sterilizer 1</b>	<b>Sterilizer 2</b>	<b>Sterilizer 3</b>	<b>Sterilizer 4</b>
<b>11-Aug</b>	Cycle 6541 Brain needles	Drill Set cable. Did not have another available. Cycle 6741	Cycle 6964 Fine cardiac instrument	Mac blade anesthesia.Cycle 6351
		Dr. Masson's Hager dilators.- one of a kind set. Cycle 6477	Cycle 6965 wire scissors	Cycle 6388 brain needles-one of a kind set
		Small bone drill immediate use. Cycle 6488	Cycle 6990 Mayo scissors	
		Cycle 6489 beaver blade handle one	Cycle 7005 Nasal fracture tray Dr. Alonso	
<b>11-Sep</b>			Cycle 7036 system cd 3	Biomet femoral peg drill
			Cycle 7037 periteneal elevator	EEA sizers- one of a kind set
				6 piece retractor hip and knee Dr. Urquhart
<b>11-Oct</b>	Cycle 7167 Lacremal Probes & DR. Morkonich instrument set.	Cycle 6601 DR. Constads chopper insts. They were not sterilized on date given	Immediate use was used no items written as to what was flashed.	Pituitary Forceps and Ronguer V Muller, Kerrisons.
	Cycle 7143 Kaiser Pilling instrument set immediate used- one of a kind set		Cycle 7123 Fine Cardiac instrument set. (one of a kind set)	Cycle 6484 synthes tap and inserter
	Cycle 7148 Fine cardiac instrument Dr. Han Han		Dr. Holmes St Marcks instrument used number of previous cases.	Cycle 6488 Bankard retractor fell to floor- not used
	Cycle 7149 right angle clamp			C=ycle 6658 Charley ret hole in wrapper evening 2nd shift
<b>11-Nov</b>		Nasal Fracture set one of a kind set	Cycle 7191Hand & ankle set (Dr. Shin's case)	Cycle 6605 light cable
		Cycle 6696Synthes canulated screw driver		Cycle 6606 bone clamp
		Cycle 6707 Dr. Alonso's dental instrument		Cycle 6588 pituitary forceps
				Cycle 6589 Metronic instrument set hole in wrap
<b>11-Dec</b>	Simms retractor Cycle 6806	Cycle 6784 I.A & hand tray (eye)	Cycle 6662 Dr. Laggars instrument set	
	Bookwalter retractor Cycle 6753	Cycle 6789 Dr. Constad eye instruments		
		Dr. Constads instruments (eye) not used cycle 6802		
<b>12-Jan</b>	Immediate use cycle was used to strilize the stryker knee scope which was damaged.	Immediate use was used to sterilize dermotome set another one of a kind set	Fine Open heart set immediate use was used for this one of a kind set.	Immedieate use was used on hip and knee instrument retractors
	Immediate use was used to also sterilize a one of a kind set named craniofix.			
<b>12-Feb</b>	Immediate use was performed to sterilize a hall micro drill cycle # 6876. this drill was used in another case and was needed again because the first set had broken and had to be placed for repair.	Immediate use was performed to on cycle 6904 to sterilze 2 dental instruments with currettes from the doctor's office.		

	Again immediate used was used to sterilize 4 loose instruments cycle #6877 and cycle #6878 for organ donation instruments			
<b>12-Mar</b>	Immediate use was used to sterilize codman neuro instruments also pituitary forceps and cloward sets # 1 and # 2 Dr. Smith	Cataract tray and I.A eye set was sterilized in immediate use cycle for eye cases.	Immediate Use was used to sterilize O/H chest retractor.	
	Cloward Instruments 2 sets.	Immediate use was used to sterilize orthogonic set for dental case.	Intestinal Clamps	
	Intestinal clamp Cycle 6967		Aortic Dilator O/H	
			Immediate used was used to sterilize O/H finochetto retractor med size was needed and there is only one	
	Biomet tray explore set had hole in wrap. Pilling chest retractor no integrator cycle 7463			
<b>12-Apr</b>	Immediate Use was used to sterilize dental extraction set	Dental extraction set again had back to back cases not enough sets available two different cases		Dental extraction set
	Dental Extraction Tray only two available. Cycle 7036	Dental extraction tray Cycle 7025		Dental extraction tray. Cycle 6893
<b>12-May</b>	None at this time	Cycle 7064 immediate use was used to sterilize the nasal fracture set	Immediate use was used but no documentation of what was placed in sterilizer was found. Only read out slip was in the envelope	
<b>12-Jun</b>	None at this time	None at this time	None	None
<b>12-Jul</b>	None at this time	Dental Extractor set	C.S. Tech no immediate use was performed. No records of Bowie test or B.I test was on record, Cs Tech was notified of this error, and was explained the importance of this test.	Immediate use was used to sterilize patella instrument set
				Immediate used was used to sterilize pituitary set.

### Jersey City Medical Center Immediate Used Audit



## Immediate Use Usage Data

	Sterilizer 1	Sterilizer 2	Sterilizer 3	Sterilizer 4	Total
<b>Aug-11</b>	1	4	4	2	11
<b>Sep-11</b>	0	0	2	4	6
<b>Oct-11</b>	4	1	3	3	11
<b>Nov-11</b>	0	3	1	4	8
<b>Dec-11</b>	2	3	1	0	6
<b>Jan-12</b>	2	1	1	1	5
<b>Feb-12</b>	2	1	0	0	3
<b>Mar-12</b>	3	2	6	0	11
<b>Apr-12</b>	2	2	0	2	6
<b>May-12</b>	0	1	1	0	2
<b>Jun-12</b>	0	0	0	0	0
<b>Jul-12</b>	0	1	0	2	3
				<b>Total:</b>	<b>72</b>